

See discussions, stats, and author profiles for this publication at: <https://www.researchgate.net/publication/273453332>

# Structural and Mechanical Characterization of Biomaterials for Tissue Engineering: Fibrous Scaffold –Quantitative Nano–Mechanical Mapping by Atomic Force Microscopy

CONFERENCE PAPER · MARCH 2015

DOI: 10.13140/RG.2.1.3765.9683

---

READS

109

5 AUTHORS, INCLUDING:



**M. Wozniak**

Warsaw University of Technology

31 PUBLICATIONS 199 CITATIONS

SEE PROFILE



**Adrian Chlanda**

Warsaw University of Technology

4 PUBLICATIONS 3 CITATIONS

SEE PROFILE



**Ewa Kijeńska**

Warsaw University of Technology

9 PUBLICATIONS 75 CITATIONS

SEE PROFILE

# Structural and Mechanical Characterization of Biomaterials for Tissue Engineering: Fibrous Scaffold - Quantitative Nano-Mechanical Mapping by Atomic Force Microscopy

Michał J. WOZNIAK<sup>\*,\*\*</sup>, Adrian CHLANDA<sup>\*</sup>, Ewa KIJENSKA<sup>\*</sup>, Wojciech SWIESZKOWSKI<sup>\*</sup> and Krzysztof J. KURZYDŁOWSKI<sup>\*</sup>

<sup>\*</sup> Faculty of Materials Science and Engineering, Warsaw University of Technology, Poland

<sup>\*\*</sup> University Research Centre – Functional Materials, Warsaw University of Technology, Poland

## Abstract

*Biodegradable polymeric fibers with nano and submicron diameters, produced by electrospinning can be used as scaffolds in tissue engineering. However is necessary to characterize their mechanical properties especially at the nanoscale. The PeakForce Quantitative NanoMechanics (PF-QNM) is recently developed AFM mode, which allows to probe mechanical properties of the material, such as: reduced Young's modulus, deformation, adhesion, and dissipation, simultaneously with topographical imaging. In this paper we are presenting results of PF-QNM characterization of two kinds of electrospun fibers: PCL and PCL/HAp. The average calculated from DMT theory Young modulus was  $3 \pm 1$  MPa for PCL mesh and  $17 \pm 3$  MPa for PCL+HAp mesh.*

**Keywords:** PeakForce Quantitative NanoMechanics (PF-QNM), Atomic force microscopy, Mechanical properties at nanoscale, Fibrous scaffold, Electrospinning.


## 1. Introduction

Biodegradable polymeric fibers, produced by electrospinning method, with nano and submicron diameters can be used as scaffolds in tissue engineering. It is possible due to their interfibrous pore size, high surface area to volume ratio, immunogenicity, biodegradability and structural similarity to the extracellular matrix (ECM) <sup>1)</sup>. In physiological conditions fibers are subjected to stresses and strains from the surrounding biological environment. Such stresses can cause permanent deformation or even failure to scaffold structure. Therefore, there is a growing need to characterize their mechanical properties, especially at the nanoscale. Atomic force microscopy (AFM) has emerged as a powerful tool in the imaging of cells and biomaterials and probing selected mechanical properties under physiological conditions <sup>2,3,4)</sup>. The PeakForce Quantitative NanoMechanics (PF-QNM) is recently developed AFM mode, which allows to measure mechanical properties of the material, such as: reduced Young's modulus, deformation, adhesion, and dissipation, simultaneously with topographical imaging <sup>5)</sup>.

## 2. Experimental Procedure

In presented work two different kinds of fibers, produced by electrospinning from solution at applied voltage of 12 kV <sup>1)</sup>. Polycaprolactone (PCL) and Polycaprolactone with 10 wt% of hydroxyapatite (PCL/HAp). The Bruker MultiMode AFM NanoScope 8 atomic force microscope and APP NANO - ACST probe were used to measure morphology, and related


nanomechanical properties of electrospun fibrous meshes. The AFM microscope were equipped with PeakForce Quantitative NanoMechanics (PF-QNM) module. The principle of PF-QNM bases on force-distance curve of interaction between AFM probe and sample <sup>2,3)</sup>, what is shown in **Figure 1**. The Derjaguin-Muller-Toporov (DMT) model were used to measure mechanical properties of the material: reduced Young's modulus, adhesion, deformation and dissipation with high spatial resolution <sup>5)</sup>. For comparizon the fibers were examined by a scanning electron microscope (Phenom Pro X, FEI SEM).



**Fig.1** Schema of force-distance (FD) curve. The curve represents interaction during approaching of AFM probe to the sample surface, contact with sample and retracting away from surface. By analyzing FD curve the mechanical properties at nano scale can be measured. Automatic analysis of curve set allows to generate maps of mechanical properties distribution simultaneously with topographical imaging.

### 3. Results and Discussion

**Figure 2** shows the micrographs of PCL and PCL/HAp electrospun fibers obtained by SEM and AFM (Height images) operating in PF-QNM mode. The DMT Modulus images show higher values of Young modulus (brighter color) for composite PCL/HAp material in comparison with pure PCL. The average values of selected nanomechanical properties are shown in the plot at the bottom of Fig. 2.



**Fig.2.** Comparison of results obtained on PCL and PCL/HAp composite. Imaging: SEM (scale bar 10 $\mu$ m), AFM, PF-QNM (DMT Young Modulus). Average results of nanomechanical properties calculated from PF-QNM images.

Mean values of Young modulus were equal:  $17 \pm 3$  MPa and  $3 \pm 1$  MPa for PCL/HAp and PCL respectively. This fact can be of course explained by ceramic addition to composite. It should be mentioned that these data were collected from polymeric meshes, so results are average for wholes constructs, and they are different than that for bulk material. However results obtained from single fibers (data not shown) confirmed these for meshes. Another important quantity, which can be calculated from PF-QNM is adhesion force between AFM probe and sample. Adhesion is extremely important for use of electrospun fibers as scaffolds for cells seeding. Cells need to adhere properly to the biomaterial to grow new tissue. It should be mentioned, that the adhesion force and other mechanical properties depends on the radius of AFM probe. In our previous work<sup>5)</sup>, we sharpened AFM probes using focused ion beam (FIB) and demonstrated that decreasing probe radius caused decreasing adhesion and increasing measured DMT modulus value. High magnification mapping of nanomechanical properties on PCL/HAp sample allowed to distinguished in easy way bioceramic filler in polymeric matrix (data not shown).

### 4. Conclusions

The PF-QNM is very promising technique for probing nanomechanical properties of biomaterials simultaneously with topographical imaging.

Mechanical properties of nanomaterials can differ from the properties of materials at macroscale, and can depend on geometry of AFM probe.

### Acknowledgements

The research were supported by grants of National Science Center of Poland under grant UMO-2011/01/B/ST8/07559

### References

- 1) E.Kijeńska, M.P.Prabhakaran, W.Swieszkowski, K.J.Kurzydowski, S.Ramakrishna, "Electrospun bio-composite P(LLA-CL)/collagen I/collagen III scaffolds for nerve tissue engineering". *Journal of Biomedical Materials Research: Part B – Applied Biomaterials* **100**, 1093-1102, (2012)
- 2) M.J. Wozniak, G. Chen, N. Kawazoe, T. Tateishi, "Monitoring of mechanical properties of serially passaged bovine articular chondrocytes by atomic force microscopy", *Micron*, **40**, 870-875, (2009)
- 3) M.J. Wozniak, N. Kawazoe, T. Tateishi, G. Chen, "Change of the Mechanical Properties of Chondrocytes during Expansion Culture", *Soft Matter*, **6**, 2462-2469, (2010)
- 4) A. Roguska, S. Hiromoto, M.J. Wozniak, M. Pisarek, A. Yamamoto, "Collagen immobilization on 316L stainless steel surface by applying electrochemical treatment: surface characterization and in-vitro study", *Applied Surface Science*, **257**, 5037–5045, (2011)
- 5) A. Chlanda, J. Rebis, E. Kijenska, M.J. Wozniak, K. Rozniatowski, W. Swieszkowski, K.J. Kurzydowski, "Quantitative imaging of electrospun fibers by PeakForce Quantitative NanoMechanics atomic force microscopy using etched scanning probes", *Micron*, *in press*